VOS LENTILLES  
MADE IN  
FRANCE

# GUIDE D'ADAPTATION AKS™

FRANCAIS

ENGLISH

ESPAÑOL

CE 1639

LCS © ALL RIGHTS RESERVED 2025 - 01/25

LCS

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
LCS

# FITTING PROCESS

- 1.** Determining the 1st intention SAG
- 2.** Lens insertion
- 3.** Clearance control after insertion
- 4.** Over-refraction
- 5.** Evaluate the 4 zones
- 6.** Stabilization axis
- 7.** Surface quality
- 8.** Lens removal
- 9.** Lens care
- 10.** Follow-up
- 11.** Managing special cases
- 12.** Modifications
- 13.** Troubleshooting
- 14.** Scleral lens decontamination instructions

# Scleral lens AKS™

Scleral lenses are large-diameter rigid lenses that extend beyond the cornea. They are designed to compensate for certain corneal pathologies, or damage to the corneal surface in cases of severe dry eye, but also in case of drop out of rigid or hybrid lenses. They can compensate for the following ametropias : Myopia, Hyperopia, Astigmatism and Presbyopia.

**MATERIAL**

Dk High (clear) or Dk200 (clear or blue)

**TREATMENT**

Tangible Hydra-PEG™ on request

**DESIGNS**

Spheric, Toric, ATD, Microvault, Multifocal and Oblate

**DIAMETER**

From 14.50mm to 18.50mm

**CENTER THICKNESS**

From 200 to 300µm for a power of -3D depending on diameter

**SAG**


From 3600 to 5600µm

**AREAS**

Cornea, Limbus and Sclera from -15 to +15 (1 step = 30µm)

**POWERS**

From -40 to +40, cylinder up to 6 and addition up to 3.50D

## 1. Determining the 1st intention SAG


Depending on the profile of the cornea and its protuberance, select the most suitable SAG.

Depending on corneal diameter and palpebral aperture, choose the most suitable diameter from the trial set.


Diameter 16.50 mm is recommended for 1st intention.

**Low elevation**  
Oblate corneas

**Normal and medium elevation**  
Surface diseases, keratoconus,  
PMD, Corneal graft  
Dry eyes

**Extreme elevation**  
Keratoglobus,  
Corneal graft

Trial set : 14 torics (diameter 16.50mm) and 10 spherics (diameter 15mm)

## 2. Lens insertion

### By the patient

Before fitting, check that the packaging and lenses are undamaged. Refer to the [instructions for use](#).

Always rinse the lens with saline or equivalent preservative-free solution, such as Cleadow SLi. Rub for a few seconds to remove any impurities and improve wettability.

#### By the practitioner

The patient tilt their head forward and pulls back the lower eyelid, while the practitioner next to him holds the upper eyelid and inserts the lens filled with preservative-free saline or Cleadow SLi (and fluorescein if required).


ECO SCLERAL

LCS SL


LCS RING


### By the practitioner


### 3. Clearance control after insertion


**NO CONTACT**  
Check the clearance


**CENTRAL CONTACT**  
Increase the SAG of +300µm

**INSUFFICIENT CLEARANCE**  
(correct at the end of the day).  
Approximately 100µm. Increase  
SAG by 200µm

**GOOD CLEARANCE AT  
INSERTION**  
Between 200 and 300µm


**EXCESSIVE CLEARANCE**  
Reduce SAG by 300 µm

### 4. Over-refraction

You can use an **autorefractometer** over the lenses, then do your conventional over-refraction, minimizing the astigmatism as much as possible. Greater than or equal to 4 dioptres, don't forget the **vertex distance conversion** when you add this correction to the lens power.

**5. Evaluate the 4 zones** (at least one hour after insertion)


**a CHECK AGAIN CENTRAL CLEARANCE**


The lens will recess by an average of 100 µm, you need a clearance slightly less than the lens thickness (between 150 and 250µm). **The slit should be as thin as possible, with high magnification and white light.**

**b CHECK THE CORNEAL AREA**


With good clearance : if you change C zone, you will increase the clearance (1 step = 30µm). The SAG parameter must be reduced accordingly.  
Example : C+3 SAG -100µm


No touch in periphery : no changes

Small peripheral touch: increase Cornea +3


Moderate peripheral touch: increase Cornea +6


Strong peripheral touch: increase Cornea +9

**c CHECK THE LIMBUS AREA**


With good clearance : if you change L zone, you increase lens clearance (1 step = 30µm). The SAG parameter must therefore be reduced accordingly.  
Example : L+3 SAG -100µm


No Limbus touch: no changes

Small touch on the Limbus: increase Limbus +3

Moderate touch on the Limbus: increase Limbus +6


Strong touch on the Limbus: increase Limbus +9

**d CHECK THE SCLERAL AREA**

Here we're going to check for compressions / blanching. These are not always visible at the insertion, and may appear as the lens recess. We recommend checking this area at the end of the day. White light + diffuse illumination.


No compression of vessels: no changes

Small compression of vessels: decrease Sclera -2


Moderate compression of vessels: decrease Sclera -4


Strong compression of vessels: decrease Sclera -6

## 6. Stabilization axis



This axis represents the « **flat** » axis of the lens, which will stabilize naturally according to the toricity of the sclera. Noted in the standard axis notation (**TABO**).

115°

15°

## 7. Surface quality

In case of poor wettability, you can request a **Tangible Hydra-PEG™\*** treatment.

Poor wettability

Good wettability

\* For the care of lenses with Tangible Hydra-PEG™ surface treatment, it is imperative to use compatible products. We therefore recommend the use of products from the LCS range or the Cleadew range, approved for Tangible Hydra-PEG™ treatment.  
Tip : for a better treatment conservation, it is possible to gently rub the lenses.

## 8. Lens removal

Use a plunger for removal. The plunger must be cleaned according to [the supplier's instructions](#).

**1**

Hold the plunger between thumb and forefinger, then look at yourself in the mirror. Position your middle finger just below the free edge of the lower eyelid and the index finger of your other hand on the free edge of the upper eyelid.

**2**

Stick the plunger on the bottom of the lens and gently pull it upwards. The lens will lose its adhesion, first from the bottom, then the whole lens will release itself. You can now pick it up with your thumb and forefinger.

HANDLING

## 9. Lens care

Always clean your lenses after wearing them, and follow the [instructions for use](#) and safety precautions.

### Trial set lenses

Use a lens surfactant to clean them for the next patient, and store them in a multipurpose rigid solution in their original case.

As a healthcare professional, you need to keep track of these lens care. LCS provides a number of documents to help you do this :

[Decontamination instructions](#)  
[Lens traceability](#)

### Wearer's lenses

Use a povidone-iodine-based solution or a peroxide hydrogen with weekly protein removal, or a multipurpose rigid solution and cleaner on a daily basis.

**LCS Clean**  
Cleaner

**LCS Perox**  
Peroxide hydrogen

**Cleodew SL**  
Povidone-iodine-based solution

**LCS Enzym**  
Protein remover tablets

**LCS GP**  
Multipurpose rigid solution

## 10. Follow-up

### 01. One week

See the wearer again with the new parameters at the end of the day, after one week of wearing the lenses. Repeat the same observations to determine whether any modifications are necessary or not.

### 02. One month

Do the same checks to validate the lenses.

### 03. One year

Check if the lenses need to be replaced, considering factors such as poor surface condition, lens deformation, or discomfort...

**Renewal**  
≤ 2 years



**Type of wear**  
Daily diurnal

## 11. Managing special cases

AKS lenses feature specific designs : oblate, quadrants at the Limbus or Sclera (ATD) and/or microvaults.

These specific features are not detailed in this fitting guide, but can be provided on request.

AKS Spheric

Tear leakage at the edge of the lens

AKS ATD

Disappearance of the leak

## 12. Modifications

**1st solution** : To have better technical advice, you can fill in this form with your observations and send it to us.

**2nd solution** : You can also use our fitting software iAdapt: AKS tab, then second tab "Lens modifications". Simply fill in the survey with your lenses observations, then click on calculate.

Annex page 24

PLEASE FILL THE FORM IN AND SEND YOUR REQUEST DIRECTLY TO YOUR DISTRIBUTOR

AKS™  
AirKone Scleral

Customer account number: \_\_\_\_\_  
Patient name: \_\_\_\_\_

OBSERVATIONS		TO INSERTION	AFTER 1HR OR MORE
RIGHT EYE	LEFT EYE		
SERIAL NUMBER OR TRIAL LENS PARAMETERS			
STABILIZATION AXIS			
CLEARANCE EVALUATION			
WRITE DOWN YOUR OBSERVATIONS			
OVER REFRACTION			

# 13. Troubleshooting


## 01.

You observe **conjunctival prolapse** (suction of the bulbar conjunctiva under the lens). You can :

- Reduce the diameter and limbal area (L-).
- Change your lens design (from toric to spheric - if the stabilization axis is horizontal).
- Do an ATD.

## 02.

If a **bubble** is present at the insertion, it will be necessary to insert again the lens (see [chapter 2](#)).

If it appears after a few minutes, this indicates a leak.

- Check the clearance, which may be excessive.
- Use fluorescein with blue light to identify the position of the leak.
- Modify lens design accordingly.


## 03.

You notice **lots of deposits** on the lens.

- Check the cleaning routine, and ensure good palpebral hygiene.
- Consider adding Tangible Hydra-PEG™ surface treatment, taking care to use only suitable products. (for example no Progent or i-clean).

## 04.

A **cloudy** clearance or **fogging** is observed. Usually, fogging is caused by a leak. We recommend wearing the lens without fluorescein. Then apply the fluorescein over the lens and observe if and where it penetrates. The next step is to resorb the leak. If no fluo penetration is observed, the liquid should be replaced with hyaluronic acid drops (as Cleadew SLi).

## 05.

The **lower part of the sclera is compressed**, and certain vessels are unable to pass under the lens. You'll need to do S-2.

## 06.

You can see that the **fluo** has only penetrated in one place, in the lower part. It will be necessary to close the S zone in this lower part, or to adjust the toricity of the lens according to the stabilization axis.

PLEASE FILL THE FORM IN AND SEND  
YOUR REQUEST DIRECTLY TO YOUR DISTRIBUTOR

Customer account number: .....

Patient name: .....

OBSERVATIONS     TO INSERTION     AFTER 1HR OR MORE

**RIGHT EYE**



**LEFT EYE**

SERIAL NUMBER OR TRIAL LENS PARAMETERS

STABILIZATION AXIS

CLEARANCE EVALUATION

WRITE DOWN YOUR OBSERVATIONS

OVER REFRACTION

## 14. Trial set scleral lens decontamination instructions

### I – Patient risk level according to instruction N° DGS/RI3/2011/449 of December 1, 2011 on updating recommendations aimed at reducing the risks of transmission of non-conventional transmissible agents during invasive procedures.

#### a) Suspect or affected patient

The diagnosis of TSE (Creutzfeldt-Jakob disease, new variant Creutzfeldt-Jakob disease, Gersann-Sträussler-Scheinker syndrome, fatal familial insomnia) must be suspected on the basis of the presence of at least one neurological clinical sign associated with intellectual or psychiatric disorders, and after ruling out any other cause. The diagnosis can only be confirmed by neuropathological examination.

These patients present the highest level of risk of iatrogenic transmission, justifying the application of maximum precautionary measures.

#### b) Patients neither clinically suspect nor with TSE

These are now grouped together :

- Patients with no particular characteristics
- Patients with individual risk factors for TSE belonging to one of the following groups :
  - History of treatment with human growth hormone
  - History of human dura mater transplantation
  - Genetic family history of a TSE linked to a mutation encoding PrP.

### II – Recommendations according to diagnosis

#### a) Suspect or affected patient

Return to LCS laboratory for destruction. Lens returns must be accompanied by the delivery note marked "return for destruction".

#### b) Patients neither clinically suspect nor with TSE

Implement the decontamination procedure mentioned in chapter III.

### III – Decontamination procedure

- Patients neither clinically suspect nor with TSE

After lens removal:

- 1) Rub the lens twice for 30 seconds with a surface-active soap.
- 2) Rinse both sides of the lens with saline or all-in-one solution for 5 seconds.
- 3) After rinsing, place the lens in its case filled with a decontamination solution for rigid contact lenses.
- 4) Store as is until next use

**Storage** : in the event of prolonged soaking, renew the soaking solution every 4 weeks.

Precautions: carefully read the precautions for use of the solutions used.

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LCS © ALL RIGHTS RESERVED 2025 - V1 01/25